

Since 1959

Asind[®] MILK-O-STRIPS[®]

FOR SURETY OF PURITY



FOR RAPID SCREENING OF
ADULTERATION
IN MILK



Parameters & LOD

Glucose	0.005%
Maltodextrin	0.05%
Hydrogen Peroxide	0.005%
Formaldehyde/	0.1%
Formalin	
Salt	0.005%
Urea	0.07%
Soda/Neutralizers	0.02%
Starch	0.01%
Soap/Detergent	0.25%
Hypochlorite	0.5%



What are Milk-O-Strip?

Milk-O-Strips are enzyme-based, ready-to-use test strips designed for the rapid detection of adulteration in milk. Whether you're a dairy farmer, processor, or quality controller, Milk-O-Strips offer a fast and accurate solution to protect purity and build consumer trust.

Why choose Milk-O-Strip?



Instant Results - No lab equipment needed



Highly Sensitive - Detects even low levels



User-Friendly - Quick dip, easy read



Trusted by Dairy Processors & Quality Teams



How it works?

Dip - Dip in the milk sample

Remove - Remove and wait

Match - Match with the provided color scale.



VALIDATED PERFORMANCE

NABL & FSSAI-notified laboratories



Sales Office: - A-101, First Floor, Group of Industries, Block - C, Wazirpur Industrial Area, Delhi - 110052, INDIA R&D Office : - Plot No - 313, Gajipur Ind. Area, Tehsil-Badhkal, Faribadab - 121001, INDIA Email Id - arflifesciencesllp@gmail.com, Contact Us - +91 - 8448294456, Website - www.arflifesciences.com